

Link do produktu: <http://www.novazym.sklep.pl/jy-zy2-western-blotting-system-trans-blot-p-360.html>



JY-ZY2 Western Blotting System (Trans-Blot)

Dostępność

Na zamówienie

Numer katalogowy

JY-ZY2

Opis produktu

JY-ZY2 Western Blotting System

Western blot system is used for testing a mixture of a great number of proteins. it can show us the information of our target protein (with comparison to a size marker or ladder in kDa), also the protein expression could be expressed by the run (with comparison to a control such as untreated sample or another cell type or tissue). "JY" series Western blot system covers almost all the areas of the general running, meeting the most people's experiments demands.

Advantages

JY-ZY2 Trans-Blot is generally intend for the transfer print, when the run finished, the Protein or Nucleic Acid will be moved to the Acetate Cellulose by using the Trans-blot system.

Key features

- Gel transfer size(W×L): 170 × 160 mm
- Buffer Volume: 4000 mL
- Recommend matched power supply: JY200C/JY600

JY-ZY2 Standard Accessory

Description	Quantity
Body cell	1
Trans-clamp	2
Gel fiber mat	4
Cooling pipe standing	1
Rubber pipe	2
Electrical Cables	1