

Link do produktu: <http://www.novazym.sklep.pl/jy-zy3-western-blotting-system-semi-dry-blotters-p-363.html>



## JY-ZY3 Western Blotting System (Semi Dry blotters)

Dostępność

**Na zamówienie**

Numer katalogowy

**JY-ZY3**

### Opis produktu

#### JY-ZY3 Western Blotting System

Western blot system is used for testing a mixture of a great number of proteins. it can show us the information of our target protein (with comparison to a size marker or ladder in kDa), also the protein expression could be expressed by the run (with comparison to a control such as untreated sample or another cell type or tissue). "JY" series Western blot system covers almost all the areas of the general running, meeting the most people's experiments demands.

#### Advantages

JY series Semi Dry blotters are designed with the selected the special processed materials, compare with general semi dry blotter, only minimum buffer needed, less heat and faster speed, compatible with the common used power supplies, safety designed, power automatically cut off by removing the lids.

#### Key features

Model	JY-ZY3	JY-ZY3+
Transfer size(mm)	140×140	240×160

Recommend matched power supply: JY300/JY300C/JY600/JY600C