

Link do produktu: <http://www.novazym.sklep.pl/jy-zy5-western-blotting-system-trans-blot-p-361.html>



## JY-ZY5 Western Blotting System (Trans-Blot)

Dostępność

**Na zamówienie**

Numer katalogowy

**JY-ZY5**

### Opis produktu

#### JY-ZY5 Western Blotting System

Western blot system is used for testing a mixture of a great number of proteins. it can show us the information of our target protein (with comparison to a size marker or ladder in kDa), also the protein expression could be expressed by the run (with comparison to a control such as untreated sample or another cell type or tissue). "JY" series Western blot system covers almost all the areas of the general running, meeting the most people's experiments demands.

#### Advantages

JY-ZY5 Trans-Blot is compatible with the JY-SCZ2, JY-SCZ2+ series vertical electrophoresis tank, meeting the general needs of protein and Nucleic Acid run, with special design and injection moulding technology, two gels can be run simultaneously under identical buffer condition, high resolution separation with simple operation. JY-ZY5 Trans-Blot adopts the safety design conception; power will be automatically cut when the lid is removed.

#### Key features

- Gel transfer size(W×L): 100 × 100 mm
- Buffer Volume: 750 mL
- Recommend matched power supply: JY200C/JY600

#### JY-ZY5 Standard Accessory

| Description       | Quantity |
|-------------------|----------|
| Body cell         | 1        |
| Trans-core        | 1        |
| Gel clamp         | 2        |
| Gel fiber mat     | 4        |
| Electrical Cables | 1        |