

Dane aktualne na dzień: 02-05-2025 01:14

Link do produktu: <http://www.novazym.sklep.pl/jy02s-uv-transilluminator-p-371.html>



JY02S UV Transilluminator

Dostępność

Na zamówienie

Numer katalogowy

JY02S

Opis produktu

JY02S UV Transilluminator

JY series UV Transilluminator is mainly applied to the electrophoresis results observation. With the principles of different wavelength shows the different electrophoresis band result, They are compatible with the observation of nuclear acid and protein electrophoresis bands, by using the high resolution of CCD and software technology, the gel pictures are clearly showed.

Specyfikaton

- Reflection ultraviolet wavelength: 254 nm 365 nm
- Source ultraviolet wavelength: 312 nm
- Transmission size: 200 × 150 mm
- Obscure camera operation could be used for 24-hour
- Equipped with gel cutting door, simple for operation
- Special camera lift built in, for adjusting the CCD position
- For Multl type gel observation