

Link do produktu: <http://www.novazym.sklep.pl/jy200c-universal-power-supply-p-197.html>



JY200C Universal Power Supply

Dostępność

Na zamówienie

Numer katalogowy

JY200C

Opis produktu

JY200C Universal Power Supply

Specification:

1. Size: 28×24×11 [cm]
2. Weight: 3.2 kg
3. Type of output: Constant-Voltage, Current or Power
4. Output Range: 5÷200 V, 1÷2000 mA, 1÷200 W
5. Increment: 1 V, 1 mA, 1 W
6. Timer Range: 1min÷99 hr, 59 min
7. Control: Intellectual microprocessor control unit

Key features:

1. Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating
2. Ten general electrophoresis methods stored function
3. Automatic memory
4. LCD screen intellectual display of Constant-Voltage, Current or Power values
5. Four sets of output jacks
6. Auto-switching function
7. Touch control panel
8. Injection molding case
9. Optional functions (Recovery of power failure, Non-proliferation electric field)