

Link do produktu: <http://www.novazym.sklep.pl/jy300-basic-power-supply-p-198.html>



## JY300 Basic Power Supply

Dostępność

**Na zamówienie**

Numer katalogowy

**JY300**

### Opis produktu

## JY300 Basic Power Supply

### **Specification:**

1. Size: 27×26×12 [cm]
2. Weight: 3.2 kg
3. Type of output: Constant-Voltage or Current (Adjustable)
4. Output Range: 5÷300 V, 1÷300 mA, 90 W (Maximum)
5. Increment: 1 V, 1 mA
6. Timer Range: 1min÷9 hr, 59 min
7. Display: Double LED digital screen display
8. Control: Intellectual microprocessor control unit

### **Key features:**

1. Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating
2. Automatic memory
3. Auto-switching function

4. LED screen intellectual display of Constant-Voltage, Current or Power values
5. Touch control panel
6. Injection molding case