

Link do produktu: <http://www.novazym.sklep.pl/jy3000-high-voltage-multi-usage-power-supply-p-203.html>



JY3000 High-Voltage Multi-Usage Power Supply

Dostępność

Na zamówienie

Numer katalogowy

JY3000

Opis produktu

JY3000 High-Voltage Multi-Usage Power Supply

Specification:

1. Size: 30×26×12 [cm]
2. Weight: 5 kg
3. Output Range: 20?3000 V, 1?400 mA, 1?400 W
4. Output jacks: Four sets of output jacks are provided to facilitate connection of up to
5. 4 identical electrophoresis cells simultaneously
6. Timer control: 1min?99 hr, 59 min
7. Display: Backlit LCD (128×64 pixel)
8. Control: Intellectual microprocessor control unit

Key features:

1. Programmable constant voltage, constant current or constant Power operation with automatic cross over
2. Safety features include automatic detection of no-load, short circuit, rapid resistance change, ground leak and system overheating

3. Ten general electrophoresis methods stored function
4. Automatic memory
5. Standard and timer operation function
6. Automatic completion of a run interrupted by power failure
7. Intellectual display of Constant-Voltage, Current, Power
8. Adjustable LCD display contrast
9. Optional functions (Recovery of power failure, Non-proliferation electric field, V-H, run step by step)