

Link do produktu: <http://www.novazym.sklep.pl/jy3000-high-voltage-multi-usage-power-supply-p-345.html>



JY3000+ High-Voltage Multi-Usage Power Supply

Dostępność

Na zamówienie

Numer katalogowy

JY3000+

Opis produktu

Specification

- Size: 37×31×13(cm)
- Weight: 9kg
- Type of Output: Constant-Voltage, Current or Power
- Output Range: 20--3000V, 1--400mA, 1--200W(double lines) Increment: 1V, 1mA, 1W.
- Timer Range: 1min--99hr, 59min
- Display: Backlit LCD(128×64pixel)
- Control: Intellectual microprocessor control unit

Key features

- Safety features include automatic detection of no-load, short circuit, rapid resistance change, ground leak and system overheating.
- Twenty general electrophoresis methods stored function.
- Gel temperature tested function (relevant temperature sensor need to be purchased separately)
- Automatic completion of a run interrupted by power failure
- Auto-switching function
- Standard and Timer operation function
- Intellectual display of Constant-Voltage, Current or Power
- Optional function: Recovery of power failure, Non-proliferation electric field, V-H?run step by step