

Link do produktu: <http://www.novazym.sklep.pl/jy600-basic-power-supply-p-200.html>



JY600 Basic Power Supply

Dostępność

Na zamówienie

Numer katalogowy

JY600

Opis produktu

JY600 Basic Power Supply

Specification:

1. Size: 27×26×12 [cm]
2. Weight: 3.2 kg
3. Type of output: Constant-Voltage or Current (Adjustable)
4. Output Range: 5÷600 V, 1÷200 mA, 120 W (Maximum)
5. Increment: 1 V, 1 mA
6. Timer Range: 1min÷9 hr, 59 min
7. Display: Double LED digital screen display
8. Control: Intellectual microprocessor control unit

Key features:

1. Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating
2. Automatic memory
3. Auto-switching function
4. LED screen intellectual display of Constant-Voltage, Current or Power values
5. Touch control panel
6. Injection molding case