

Link do produktu: <http://www.novazym.sklep.pl/lumi-single-tube-luminometer-p-882.html>



Lumi - Single Tube Luminometer

Opis produktu

Lumi - Single Tube Luminometer



Lumi is a Single Tube Luminometer featuring dominant sensitivity through PMT Technology. It is not only convenient featuring simple screen, but also easy to use due to user-friendly User Interface with various detection modes and setups.

design with 7.0" LCD touch

?Single Tube Measurement

?High Sensitivity by PMT Technology

?Compact Design

?7.0" LCD Touch Screen

?Easy Data Saving by USB

?Various Setup Variables

```

graph LR
    A[Power on] --> B[Insert the sample tube]
    B --> C[Set up the Modes and Variables]
    C --> D[Start measurement]
    D --> E[Confirm the results]
    E --> F[ ]
  
```

The flowchart illustrates the sequential steps for operating the HemoCue 4610 device. It begins with 'Power on', followed by 'Insert the sample tube', 'Set up the Modes and Variables', 'Start measurement', and 'Confirm the results'. The final step is an empty box, indicating the end of the process or a transition to the next phase.

BRANK
ZDJĘCIA

