


Microtube centrifuge HWLAB

Opis produktu



Description

HWLAB® Microtube centrifuge is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Product features

- Maintenance-free brushless motor
- LCD display of time and speed
- Adjustable speed from 500r/min to 8500r/min
- Timer, from 15s to 99min, or continuous operation
- 12×0.5ml + 12×1.5ml/2.0ml + 4×8×0.2ml capacity rotor
- Lid automatically open at the end of the run to prevent sample warming and to reduce handling time

Technical specifications

Description	HW56
Maximum rcf	4850xg
Power supply	AC100V-240V/50-60Hz
Rotor capacity	12 x 0.5ml + 12 x 1.5ml/2.0ml + 4 x 8 x 0.2ml
Speed	From 500r/min to 8500r/min
Timer	15s-99min / continuous operation
Motor type	Brushless motor
Dimension	212x270x158mm
Weight	4.3kg



Ordering information: HWLAB® Microtube centrifuge

	Model	Cat.NO	Destription
	HW56	C-10000	Widely ues in clinical labs and donation centers
	HW56.1	C-10001	For 12x0. 5ml, 12x1. 5ml/2. 0ml,4x8x0.2ml tubes