


Link do produktu: <http://www.novazym.sklep.pl/mini-microplate-centrifuge-hw4-hwlab-p-3813.html>



Mini Microplate Centrifuge HW4 HWLAB

Opis produktu



Description

HWLAB[®] mini microplate centrifuge is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Technical specifications

Description	HW4
Maximum rcf	580xg
Power supply	AC100V-240V/50-60Hz
Rotor capacity	8x96-well plates / 4xdeepwell plates
Speed	From 500r/min to 2500r/min
Timer	15s-99min / continuous operation
Motor type	Brushless motor
Dimension	285x225x190mm
Weight	4kg

Product features

- Maintenance-free brushless motor
- LCD display of time and speed
- Adjustable speed from 500r/min to 2500r/min
- Timer, from 15s to 99min, or continuous operation
- Wide applicable range. Suitable for skirted and non-skirted of all standard of PCR plates and maximum loading of 28mm deepwell plates
- Super capacity, it can running for 8 PCR 96-well plates or 4 deepwell plates
- Humanized design. Three operating modes can choose, including short spin/ setting time/continuous operation
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time



Ordering information: HWLAB[®] Mini microplate centrifuge

Model	Cat.NO	Destription
 HW4	C-15000	Widely ues in clinical labs and donation centers