

Link do produktu: <http://www.novazym.sklep.pl/mini-vortex-mixer-p-3657.html>



Mini Vortex Mixer

Opis produktu

Mini Vortex Mixer



Features:

- ? Small footprint and compact structure
- ? Touch function
- ? Fixed speed of 3000rpm
- ? Orbital diameter of 5.5mm
- ? Patented design of foot pad can hold the vortex mixer firmly in place when in mixing
- ? Full containers of 50ml can be activated for mixing within 5 seconds
- ? Power supply has been CE, cULus, CCC, SAA, FCC and GS certified

Specifications:

Model	MI0101001
Rated voltage [V AC]	100-240
Rated Frequency [Hz]	50-60
Output Voltage [V DC]	12
Output Power [W]	7.2
Speed [1/min]	3000
Operating time [%]	100
Protection Class	II
Pollution Degree	2
Overvoltage Category	II
EN 60529 Protection Class	IP 42
Ambient Humidity (Relative) [%]	80
Ambient Temperature [°C]	5 to 40 (in the operation)
Transport Temperature [°C]	-25 to +70
Storage Temperature [°C]	-10 to +70
Weight [kg]	0.55
Test Tube Size [mm]	Φ30
Dimensions: (Φ x H) [mm]	100 x 78
Type of Movement	Orbital
Orbital Diameter [mm]	6
Touch Function	Yes