

Link do produktu: <http://www.novazym.sklep.pl/miniwirowka-laboratoryjna-15k-o-wysokiej-predkosci-high-speed-mini-centrifuge-p-256.html>



Miniwirówka laboratoryjna 15K o wysokiej prędkości (High-speed Mini Centrifuge)

Dostępność

Na zamówienie

Numer katalogowy

MC-15K

Opis produktu

High-speed Mini Centrifuge 15K

Mini-15K high speed centrifuge is a powerful and easy-to-use benchtop assistant, under such compact design, the max. speed could be up to 14,500 rpm, and the premium brushless motor enables extremely quiet operation



Features

- ▶ Extremely quiet operation with LCD display
- ▶ Capacity: 12 x 1.5/2.0ml tubes.
- ▶ Separate short-spin key
- ▶ Low sample heating temperature (only 10°C after 30mins, at max. Speed)
- ▶ Autoclavable rotor (121°C, 20mins)
- ▶ Fast accelerate and decelerate speed (within 15 seconds)
- ▶ Metal made rotor, inner centrifuge lid and lid latch to ensure extra safe operation
- ▶ Automatic rpm/rcf conversion

Rotor and other accessories

1. High intensive clear PC cover
2. Autoclavable rotor (121°C, 20 mins)
3. Optional adapters
 - Adapters for 0.2ml tubes
 - Adapters for 0.5ml tubes



Specifications

Model	Mini-15K
Power	100 W
Max. rotational speed	14,500 rpm
Relative centrifugal force (RCF)	14,000 × g
Timer	20 sec ~ 99 min
Acceleration time to max speed	15 sec
Braking time from max speed	15 sec
Dimensions (D x W x H) cm	27 × 23 × 15.5 cm
Weight (with rotor)	4.8kg

Order informations

Code	Description
AS-08110-00	Mini-15K AC100V or 240V 50/60 Hz
AS-08041-01	Adapters for 0.5ml tubes
AS-08041-02	Adapters for 0.2ml tubes