

Link do produktu: <http://www.novazym.sklep.pl/miniwirowka-laboratoryjna-6kc-mini-centrifuge-p-458.html>



Miniwirówka laboratoryjna 6KC (Mini-Centrifuge)

Dostępność

Dostępny

Numer katalogowy

MC-6KC

Opis produktu

Mini-6KC Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a "personal" mini centrifuge. It includes 2 quick release interchangeable rotors, two speed options (4,000rpm & 6000rpm). The electronic brake can reach quick deceleration to reduce handling time.

Features

1. New rotor design for increased sample capacity
2. Quick release rotor hub allows for easy rotor exchange
3. The electronic braking system for rapid sample processing
4. Two speed options: 4,000 & 6,000rpm

Specification

1. Relative centrifugal force: 850 x g (4,000rpm) / 2,000 x g (6,000rpm)
2. Sample capacity: 8x1.5/2.0 ml tube rotor & 4 place PCR Strip rotor
3. Speed of rotation: 4,000rpm & 6,000rpm
4. Maximum noise: ≤56Db
5. Dimensions: 135 x 147 x 107mm (D x W x H)

Two speed options:

BRAK
ZDJĘCIA



Two rotors:

BRAK
ZDJĘCIA



Code	Description
AS-08020-00	Mini-6KC AC100V~240V 50~60Hz 1A
AS-08021-01	8-1.5/2.0ml angle rotor
AS-08021-02	4 -8x0.2ml strip rotor
AS-08061-05	Adapters for 0.5~0.6ml tubes
AS-08061-06	Adapters for 0.2ml tubes