

Link do produktu: <http://www.novazym.sklep.pl/miniwirowka-laboratoryjna-6ks-mini-centrifuge-p-459.html>



Miniwirówka laboratoryjna 6KS (Mini-Centrifuge)

Dostępność

Dostępny

Numer katalogowy

MC-6KS

Opis produktu

Mini-6KS Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a "personal" mini centrifuge. It is supplied with one rotor suitable for 6-1.5/2.0ml tubes, 6-0.5ml tubes and 2-8x0.2ml PCR tubes, and two speed options (4,000 & 6,000rpm). The electronic braking can reach quick deceleration to reduce handling time.

Features

1. Patented 3-in-1 rotor design, high capacity and time-saving
2. Rotor -- 6-1.5/2.0ml, 6-0.5ml tubes and 2-8x0.2ml PCR strip tubes
3. The electronic braking system for rapid sample processing
4. Two speed options: 4,000 & 6,000rpm

Specification

1. Relative centrifugal force: 750 x g (4000rpm) / 1700 x g (6000rpm)
2. Sample capacity: 6x1.5/2.0ml, 6x0.5ml tubes and 2-8x0.2ml PCR Strip rotor
3. Speed of rotation: 4000rpm & 6000rpm
4. Maximum noise: ≤56Db
5. Dimensions: 135 x 147 x 120mm (D x W x H)

Two speed options:



Code	Description
AS-08030-00	Mini-6KS AC100V~240V 50~60Hz 1A
AS-08031-01	6x1.5/2.0ml, 6x0.5ml tubes and 2-8x0.2ml PCR Strip rotor