

Link do produktu: <http://www.novazym.sklep.pl/miniwirowka-laboratoryjna-6kt-mini-centrifuge-p-460.html>



Miniwirówka laboratoryjna 6KT (Mini-Centrifuge)

Dostępność

Dostępny

Numer katalogowy

MC-6KT

Opis produktu

Mini-6KT Mini-Centrifuge can be supplied with different rotors for different applications. It is also offered two speed options (4,000 & 6000rpm). The electronic braking can reach quick deceleration to reduce handling time.


Features

1. Patented rotor design for 48well 50ul PCR plate
2. Optional rotor: 4 x 5ml Eppendorf-tube rotor
3. The electronic braking system for rapid sample processing
4. Two speed options: 4,000 & 6,000rpm

Specification

1. Relative centrifugal force: 850 x g (4,000rpm) / 2,000 x g (6,000rpm)
2. Sample capacity: 48well 50ul PCR plate rotor
or 4 x 5ml Eppendorf-tubes rotor
3. Speed of rotation: 4,000rpm & 6,000rpm
4. Maximum noise: ≤56Db
5. Dimensions: 135 x 147 x 120mm (D x W x H)

Two speed options:



Rotors:



48well 50ul PCR plate rotor



4 x 5ml tubes rotor (Optional)

Code	Description
AS-08100-00	Mini-6KT AC100V~240V 50~60Hz 1A
AS-08101-01	48well 50ul PCR plate rotor
AS-08101-02	4 x 5ml tubes rotor (Optional)
AS-08101-03	Adaptors for 5ml normal tubes (Optional)