

Link do produktu: <http://www.novazym.sklep.pl/miniwirowka-laboratoryjna-7k-mini-centrifuge-p-253.html>



Miniwirówka laboratoryjna 7K (Mini Centrifuge)

Dostępność

Na zamówienie

Numer katalogowy

MC-7K

Opis produktu

[Mini-6K](#) / [Mini-7K](#) / [Mini-10K Centrifuge](#)

Model	Mini-6K	Mini-7K	Mini-10K
Power requirement	30 W	30 W	40 W
Max. rotational speed	6 000 rpm	7 000 rpm	10 000 rpm
Relative centrifugal force (RCF)	2 000 x g	2 300 x g	5 000 x g
Rotors	6 x 1.5/2.0 mL angle rotor	6 x 1.5/2.0 mL angle rotor	6 x 1.5/2.0 mL angle rotor
	2 x 8 x 0.2 mL strip rotor	2 x 8 x 0.2 mL strip rotor	
	Rotor for slide		

BRAK
ZDJĘCIA



BRAK
ZDJĘCIA



BRAK
ZDJĘCIA

