

Link do produktu: <http://www.novazym.sklep.pl/miniwirowka-laboratoryjna-plytkowa-mini-p25-mini-plate-centrifuge-p-457.html>



Miniwirówka laboratoryjna płytkowa - Mini P25 (Mini- Plate Centrifuge)

Dostępność

Dostępny

Numer katalogowy

Mini-P25

Opis produktu

The compact and economical Mini-P25 Mini-Plate Centrifuge with two-position rotor is available from pacific laboratory products. It can quickly spin down droplets and condensation and is for use before and after thermal cycling to improve PCR yield. It accepts skirted, non-skirted and all standard PCR plates making, and it is suitable for the molecular laboratory.

Features

1. Quickly spin down droplets
2. Use before and after thermal cycling to increase PCR yield
3. Accepts skirted, non-skirted and all standard PCR plates
4. Less than 1/4 the size of most plate centrifuges

Specification

1. Capacity: 2 PCR Plates 96tw
2. Rotor Speed: 2500rpm
3. G Force: 500×g
4. Dimensions: 20(W)×23(D)×19.3(H)cm
5. Weight: 3kg
6. Power: 120V/230V, 50-60HZ

Operation

1. Open the lid and adjust the rotor
2. Put the PCR plates in the rotor
care about the direction

BRAK
ZDJĘCIA



BRAK
ZDJĘCIA



3. Close the lid. Switch on the button
hold the button for 10-15 seconds

BRAK
ZDJĘCIA



4. Then switch off the button.

BRAK
ZDJĘCIA

