

Link do produktu: <http://www.novazym.sklep.pl/miniwirowka-laboratoryjna-plytkowa-mpc-22-mini-plate-centrifuge-p-456.html>



Miniwirówka laboratoryjna płytkowa MPC-22 (Mini-Plate Centrifuge)

Dostępność	Dostępny
Numer katalogowy	MPC-22

Opis produktu

Centrifuge MPC-22

MPC-22 Mini-Plate Centrifuge is the first one adopted an unique design for the horizontal rotor centrifuge (Patented) in the world. The plate or tubes set into the rotor with 2pcs slope-rotor-shelf. It can quickly spin down droplets and condensation and is for use before and after thermal cycling to improve PCR yield and other tests. The MPC-22 combines functionality, easily operation and comfort ? a lot of ?high? qualities to come standard in one compact micro-centrifuge

Features

1. Unique design for the horizontal rotor (Patented)
2. It can accept below formats:
 - ? PCR standard plates with skirted, non-skirted
 - ? Elisa Plates
 - ? 8x0.2ml PCR strip tubes
 - ? 0.2ml PCR tubes
3. Braking time from max.speed: lower than 30sec
4. Timer can be set up to 10minutes
5. Short-spin and timing of two operations make the tests more convenient
6. Streamline design with more compact shape and transparent PC lid

Specification

1. Capacity: 2pcs PCR plates with skirt or without skirt / Elisa plate, 8x0.2ml PCR strip tubes / 0.2ml PCR tubes
2. Rotor Speed: 2,200rpm
3. G Force: 480×g
4. Timer: 0~10min
5. Dimensions: 290 x 360 x 140mm (W x D x H)
6. Weight: 3.6kg
7. Power: AC120V/220V, 50-60Hz

Functional diagram

BRAK
ZDJĘCIA



Before working

working

Adaptors for MPC-22

BRAK
ZDJĘCIA



Adaptor for PCR plates with
skirt and without skirt

BRAK
ZDJĘCIA



Adaptor for 8x0.2ml strip or
0.2ml PCR tubes

Code	Description
AS-08090-00	MPC-22 AC120V or AC220V 50~60Hz 1A
AS-08091-01	Adaptor for PCR plates with skirt and without skirt
AS-08091-02	Adaptor for 8x0.2ml strip or 0.2ml PCR tubes