

Dane aktualne na dzień: 17-05-2024 23:53

Link do produktu: <http://www.novazym.sklep.pl/multi-function-mini-centrifuge-hwlab-p-3807.html>



Multi-function mini centrifuge HWLAB

Opis produktu

Multi-function mini centrifuge HWLAB



Description

HWLAB® Multi-function mini centrifuge is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Technical specifications

Description	HW6
Maximum rcf	6900xg
Power supply	AC100V-240V/50-60Hz
Rotor capacity	6 x 1.5ml/2.0ml, 2 x 8 x 0.2ml, 4 x 5ml, 12 x 30mm/40mm
Adapter	6 x 0.2ml/0.5ml, 4 x 1.3ml
Speed	From 1000r/min to 12000r/min
Timer	15s-99min / continuous operation
Motor type	Brushless motor
Dimension	200 X 160 X 130mm
Weight	0.9kg

Product features


- Maintenance-free brushless motor
- LCD display of time and speed
- Adjustable speed from 1000r/min to 12000r/min
- Timer, from 15s to 99min, or continuous operation



Rotor #1



Rotor #2



Rotor #3




Rotor #4



Rotor #5

Ordering information: HWLAB® Multi-function mini centrifuge

Model	Cat.NO	Destription
 + 	L-2109	HW6 centrifuge with rotor#1 and rotor#2 for 6x1.5ml/2.0ml tubes and 2x8x0.2 PCR tubes, adapter for 0.2ml/0.5ml tubes.
 + 	L-2109A	HW6 centrifuge with rotor#3 for 4x5ml tubes, adapter for 1.3ml tubes.
 + 	L-2109B	HW6 centrifuge with rotor#4 for 12x30mm/40mm hematocrit tubes.
 + 	L-2109C	HW6 centrifuge with rotor#5 for 6x1.5ml/2.0ml tubes and 6x0.5ml tubes as well as 2x8x0.2 PCR tubes.
 40mm hematocrit tube	L-2109B1	Suitable for rotor#4, 40mm length, plain, non-heparinized, 200pcs/vial, 20 vial/box.