


Link do produktu: <http://www.novazym.sklep.pl/multi-function-shaker-mx1-hwlab-p-3815.html>



Multi-function Shaker MX1 HWLAB

Opis produktu



Description

HWLAB® Multi-function Shaker suitable for shaking tasks with all small vessels and microplates, widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Product features

- Suitable for shaking tasks with all small vessels and microplates
- LCD display of time and speed
- Wide range of attachments
- Timer with countdown function
- Adjustable speed from 300r/min to 3000r/min
- Timer, from 1s to 24h or continuous operation
- Wide voltage, wide frequency, AC 100V-240V, 50-60Hz
- Maintenance-free brushless motor
- Shaker diameter 4.5mm

Technical specifications

Description	MX1
Voltage	AC 100-240V
Frequency	50-60Hz
Speed range	300r/min-3000r/min
Timer	1s-24h / continuous operation
Shaking motion	Orbital
Orbital diameter	4.5mm
Display	LCD display speed and time
Motor type	Brushless motor
Dimension	210 X 150 X 65mm
Weight	2.5kg

Ordering information: HWLAB® Multi-function Shaker

Model	Ord.No	Description
 MX 1 Multi-function Shaker	H-10000	Suitable for mixing tasks with all small vessels and microplates
 MX 1.1 Standard attachment	H-10001	For test tubes and small vessels upto Ø 50 mm
 MX 1.2 Round attachment	H-10002	For different rubber foam inserts
 MX 1.3 Microplate attachment	H-10003	For one microplate
 MX 1.2.1 Microdish attachment	H-10004	For inserting into the round attachment
 MX 1.2.2 Test tube insert	H-10005	You can make holes as you like
 MX 1.2.3 Test tube insert	H-10006	For 14 test tubes Ø 10 mm
 MX 1.2.4 Test tube insert	H-10007	For 6 test tubes Ø 12 mm
 MX 1.2.5 Test tube insert	H-10008	For 4 test tubes Ø 16 mm