

Link do produktu: <http://www.novazym.sklep.pl/multi-tube-vortexer-model-mtv-100-p-248.html>



Multi-tube Vortexer, Model MTV-100

Dostępność

Na zamówienie

Opis produktu

Specification

The MTV-100 Multi-Tube Vortex Mixer, mix one can handle up to 50 samples. Speed is adjustable from 500 to 2500 rpm with electronic speed control. Built-in timer is adjustable from one minute to 99 h 59 min operation. Various tube sizes can be mixed by simply changing the foam tube racks.

Programmable pulsing mode allows user to adjust the pulse-on and pulse-off times between 1 ~ 59 seconds in 1 second intervals. This feature enhances the vortex action by creating a more vigorous mix.

Speed range: 500 - 2500 rpm

Speed accuracy: ± 50 rpm (at 2 500 rpm detection)

Timing range: 1 min-100 h

Working diameter: 3.6 mm

Maximum load: 4.5 kg

Dimensions: 250 x 426 x 480 mm

Weight: 21.5 kg

Power: 230 V or 120 V, 50/60 Hz, 60 W

Key features

- Can handle 50 samples
- Two modes: continuous type of work or Programmable pulsing mode
- Microprocessor control
- LCD display, and easy to use the panel
- The variable speed microprocessor control provides consistent uniform mixing action
- Independent LCD displays for speed and time allow operator to view both settings at once