

Link do produktu: <http://www.novazym.sklep.pl/nabi-spektrofotometr-p-4132.html>



Nabi - spektrofotometr

Opis produktu

Nabi - spectrophotometer

MicroDigital's Nabi, UV/Vis Nano Spectrophotometer, measures both cuvette type and micro-volume sample with supreme accuracy, outstanding precision, and superb repeatability.

Applications:

- DNA/RNA measurement: dsDNA, ssDNA, RNA, and Oligonucleotide.
- Measurement on Protein samples: Bradford, Lowry, BCA, Biuret, DirectUV and user's own wavelength setup.
- Measurement on Cell Culture samples.
- Kinetic Measurement.
- End Point Measurement.

Parameters:

Project	Parameters
Range	200 - 1,100nm
Light source	Xenon flash lamp
Detector	CCD (2,048 Pixels)
Minimum Sample Volume	0.5 μ l
Minimum Measurement Time	3 Seconds
Absorbance Accuracy	$\pm$ 0.002 abs.
Spectral Resolution	0.3 nm
Detection Limit	4 ng / μ l(dsDNA)
Maximum Conc.	10 μ g / μ l(dsDNA)
LCD Touch Panel	7.0 inch
Dimension	W : 220mm, D : 280mm, H : 220mm
Power	100-240V~, 50/60Hz
Weight	4.8 kg
Warranty	1 Year


Packing:




- Main instrument
- Pedestal cover
- Cuvette stick
- Cuvettes
- Power cord
- User manual

BRAK
ZDJĘCIA

Main instrument:





Workflow: