

Link do produktu: <http://www.novazym.sklep.pl/agile-plate-illuminator-p-320.html>



Agile Plate Illuminator

Dostępność	Na zamówienie
Czas wysyłki	14 dni
Numer katalogowy	APL-96

Opis produktu

Agile Plate Illuminator

Product specification:

Agile Illuminator provides a simple way to track the loading of reactants by illuminating the target wells which can be single well or a complete/partial row or column of wells, or all 96 wells.

Key features:

- Compatible for 96-well PCR Plate and ELISA Plate
- 3 selectable light intensities
- Applicable with all standard single channel, 8-channel, 16-channel pipettors
- Intuitive and easy instrument with pre-programmed loading patterns
- Water resistant, Compact and Portable

Art. Nr

APL-96

Size

18.8 (L) x 14 (W) x 1.2 (H) cm

Weight

280 g

Vessel Types

96-well PCR Plate and ELISA Plate

Programs

12 preset and user-selectable loading patterns

Light Source

White LED or Red LED (optional)

Operating Voltage

3.5 . 5.3V

Operating Power Consumption

400 mA

Power Supply

100 - 240 V input power, 5 VDC, 1A

Package

Agile Illuminator Unit, Foot pedal and mini-USB