


Link do produktu: <http://www.novazym.sklep.pl/bandpeeper-for-mupid-exu-series-real-time-electrophoresis-p-391.html>



BandPeeper for Mupid exU series (Real-Time Electrophoresis)

Dostępność

Na zamówienie

Numer katalogowy

RT-172

Opis produktu

[BandPeeper for Mupid exU series \(Real-Time Electrophoresis\) MaestroGen](#)

MaestroGen offers BandPeeper real-time electrophoresis LED illuminator for DNA gel fluorescent signal observation in water solution.

[Specification](#)

Size: 166*170*51 mm
Window size: 150*60 mm
Filter type: Amber filter
Filter surface: Anti-Fog treatment
PCBA surface: Anti-moisture treatment
Power: 24V - 1.25A
Sensitivity: ~ 2 (ng)
LED wavelength: 470 nm
Compatibility: only with Mupid tanks (Mupid exU series: exU, ACE, and ONE)
Weight: 0.9 kg

[Key features](#)

- Time saving for new sample electrophoresis trial & error
- Cost saving for repeated trial & error
- Accurate electrophoresis duration control
- Compatible with Mupid-size tank (Mupid exU series: Mupid exU, Mupid ACE, and Mupid ONE)
- Strong blue LED: used with common nucleic acid dyes
- High sensitivity (~2 ng)
- Anti-fog design on amber filter; anti-water condensation
- No decay for DNA gel (UV flashlamp decays DNA gel)
- Competitive price