

Link do produktu: <http://www.novazym.sklep.pl/microplate-incubator-mk100-4p-p-308.html>



Microplate Incubator MK100-4P

Dostępność

Na zamówienie

Numer katalogowy

MK100-4P

Opis produktu

Microplate Incubator MK100-4P

Specification

The Incubator MK100-4P is a high performance microplate incubator which accommodates up to four microplates. The MK100-4P can be used for any enzyme or cell-based assays requiring uniform and strictly controlled incubation up to 60 °C.

Incubation	Incubation range from 8 °C above ambient to 60 °C
	0.1 °C steps
	Temperature stability ± 0.3 °C
	Temperature gradient < 0.5 °C across the plate (@37 °C)
Warming Up Speed	< 25 min from 20 °C to 60 °C
Capacity	Up to 4 microplates
	Max. microplate dimensions: Width: 8.9 cm, Length: 13.3 cm, Height: 2.3 cm
Power Requirements	115/230 V, 60/50 Hz, max. 90 VA

Dimensions	Width: 38.5 cm, depth: 39.5 cm, height: 19 cm
	Weight: 7 kg

Key features

- Microprocessor controlled time and temperature
- Beep-signal and stop of after program completion
- Heating plates above and below the microplates
- Excellent temperature accuracy and uniformity across the plate
- User definable parameters - temperature and time
- Platform capacity: 4 microtest plates