

Dane aktualne na dzień: 21-09-2024 00:30

Link do produktu: <http://www.novazym.sklep.pl/centrifuge-for-microplates-mpc2800-p-4011.html>



Centrifuge for Microplates MPC2800

Opis produktu

Centrifuge For Microplate

MPC2800



Description:

- ◆ MPC2800 is designed for quick spins of samples in PCR plates. Primarily designed for use with 96 or 384 well plates and small volumes, it will accommodate skirted, non skirted and semi skirted styles.

Features:

- ◆ Easy to use, real-time display of all operating information.
- ◆ Countdown mode support
- ◆ DC brushless motor, low noise, small interference, maintenance-free.
- ◆ When opening the cap will automatically stop running.
- ◆ Fast acceleration from standstill to maximum speed fully in just 10 sec.
- ◆ Braking time is short, from the highest speed to a complete standstill in just 4 sec.

Parameters:

Project	Parameters	Project	Parameters
Rotating speed	2200~2800Rpm	Lid Open Auto Stop	Yes
Timing range	15sec ~1min59sec	Short Run	Yes
Max(RCF)	550g	Capacity	2 PCR plates
Acceleration time	≤10sec(0 to 2800Rpm)	Power	45W
Brake time	≤4sec(2800 to 0Rpm)	Voltage	AC110~220V/50-60Hz
Dimensions	260x240x226mm	Fuses	250V/3A Φ5×20
Weight	3.3Kg	---	---