

Link do produktu: <http://www.novazym.sklep.pl/przeciwciala-poliklonalne-polyclonal-antibody-service-pakiet-i-p-439.html>



# Przeciwciała poliklonalne (polyclonal antibody service) - PAKIET I

Dostępność

**Dostępny**

Numer katalogowy

**PL001-01**

## Opis produktu

### Custom Polyclonal Antibody Production:

This package is designed for generating antibodies against a client-supplied antigen. The package uses a client-provided antigen, and includes an 8-week immunization protocol in a species of your choice. We can perform services (ELISAs, purifications, western blotting, IHC etc.) for any project to assist you in evaluating the antibodies produced.

|                                   | Delivery                                            | Total Time  |
|-----------------------------------|-----------------------------------------------------|-------------|
| Client-provided antigen           | 2-5mg lyophilized affinity-purified antibodies      | 15-16 weeks |
| 2 rabbits or 1 goat               | Lyophilized pre-immune serum                        |             |
| Pre-Immune Serum from each animal | Elisa reports of 3 <sup>rd</sup> immunization serum |             |
| 4 immunizations per animal        | Elisa reports of post-affinity antibodies           |             |
| Production Bleeds                 |                                                     |             |
| Peptide affinity purification     |                                                     |             |
| Elisa test                        |                                                     |             |