

Dane aktualne na dzień: 01-05-2025 21:08

Link do produktu: <http://www.novazym.sklep.pl/przenosne-molekularne-laboratorium-diagnostyczne-p-3975.html>



Przenośne Molekularne Laboratorium Diagnostyczne

Numer katalogowy

MDB-1

Opis produktu

Novazym Laboratory Instruments

The Portable Molecular Diagnostics Lab



The MD-Box-Lab

The MD-Box-Lab is the only nucleic acid analysis case that can be claimed as regular bag on flights. Moreover, all devices inside the box weigh 1-2 kg and therefore are portable. Furthermore, the devices are compatible with a battery pack or car charger. The MD-Box-Lab provides all you need for molecular diagnostic analyses including devices, consumables, and protective suits.

Application:

1. Clinical diagnostics
2. Animal disease detection
3. Food safety
4. Scientific research

Features:

1. Complete solution for mobile molecular diagnostics – devices, consumables, and protective suit.
2. Portable – 12 V DC, compatible with battery pack and car charger.
3. User friendly – easy handling, quick and reliable results.
4. Flexible solution – fill-up with your different needs.

Content of MD-Box-Lab:

Description	Model	Amount
<i>Tissue Grinder</i>	G10	1
<i>Pipette</i>	0.1-2.5 µl	1
<i>Pipette</i>	2-20 µl	1
<i>Pipette</i>	20-200 µl	1
<i>Pipette</i>	100-1000 µl	1
<i>Magnetic Separator</i>	8x 1.5/2.0 ml	1
<i>Mini Centrifuge</i>	C1008	1
<i>Real-Time PCR Cyclor</i>	Mini-8	1
<i>Battery Pack</i>	Batt-3	1
<i>Dry Bath (with battery)</i>	H2O3-100CB	1
<i>Block</i>	QBLOCK	4
<i>PC Tablet</i>	Lenovo YOGA Tablet 2-1051F	1
<i>User Manual</i>	–	1



