

Link do produktu: <http://www.novazym.sklep.pl/single-speed-mini-centrifuge-hwlab-p-3804.html>



Single speed mini centrifuge HWLAB

Opis produktu

Single speed mini centrifuge HWLAB



Description

HWLAB® single speed mini centrifuge is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Product features

- Unique shape design
- Unique open button
- Low noisy
- Wide voltage, wide frequency

Technical specifications

Description	4k	6k	10k
Maximum rcf	770xg	1700xg	4800xg
Speed	4000r/min	6000r/min	10000r/min
Rotor capacity	6 x 1.5ml/2.0ml, 2 x 8 x 0.2ml		
Adapter	6 x 0.2ml/0.5ml		
Power supply	AC100V-240V/50-60Hz		
Dimension	17 x 15 x 13cm		
Quantity per carton	4pcs		
G.W.	6.0kg		
Dimension per carton	44 x 40 x 36cm		



Rotor #1



Rotor #2

Ordering information: HWLAB® Single speed mini centrifuge

Cat.NO	Description
L-2096	4k mini centrifuge, incl. Rotor #1 and Rotor #2, adapter for 0.2ml/0.5ml tubes
L-2097	6k mini centrifuge, incl. Rotor #1 and Rotor #2, adapter for 0.2ml/0.5ml tubes
L-2098	10k mini centrifuge, incl. Rotor #1 and Rotor #2, adapter for 0.2ml/0.5ml tubes
L-2102	Rotor #1 for 1.5ml/2.0ml tubes
L-2103	Rotor #2 for 2 x 8 x 0.2ml PCR tubes
L-2107	Adapter for 0.2ml/0.5ml tubes