

Link do produktu: <http://www.novazym.sklep.pl/slimwhite-transilluminator-white-light-panel-easy-to-view-colorimetric-gels-p-382.html>



SlimWhite Transilluminator (white light panel, easy to view colorimetric gels)

Dostępność

Na zamówienie

Numer katalogowy

SW-01

Opis produktu

[SlimWhite Transilluminator \(SW-01\)](#)
[MaestroGen](#)

[Specification](#)

Model Name	SW?01
Wavelength (nm)	White (blue chip + yellow fluorescent powders)
Light source	LEDs
Dimensions (mm)	210 * 305 * 15
Illumination area (mm)	170 * 275
Battery box	Optional
Weight (kg)	1.3
Power	12 VDC, 0.5 A

[Key features](#)

- Uniform illumination design (lightguide plate + diffuser)
- Cost Effective
- For viewing protein gels, coomassie blue or silver stained media
- Compact size for lab
- Low power consumption
- No need to use UV for white light conversion
- Long lifetime (~50,000 hrs)
- CE/FCC Certified
- [white light panel - easy to view colorimetric gels !!!](#)