

Link do produktu: <http://www.novazym.sklep.pl/system-do-dokumentacji-zeli-z-transiluminatorem-led-slide-imager-with-led-transilluminator-p-390.html>



System do dokumentacji żeli z transiluminatorem LED (Slide Imager with LED Transilluminator)

Dostępność

Na zamówienie

Numer katalogowy

SML-01

Opis produktu

Normal 0 21 false false PL X-NONE X-NONE MicrosoftInternetExplorer4 [Slider Imager with LED](#)

[Transilluminator MaestroGen](#)

MaestroGen provides Slider with LED transilluminator for gel photography imaging and cutting.

[Key features](#)

- Easy: gel photography imager and cutting at the same instrument
- Affordable: competitive price
- Quality: branded digital camera, ~100 % dark room for photography
- Friendly: USB interface, no need to take out and insert memory card frequently
- Considerate: rotatable LCD monitor for user to easily perform photography
- Safe: cutting gel with MaestroGen 470 nm LED transilluminator
- Easy Powering: NO need to change and re-charge batteries on camera