

Link do produktu: <http://www.novazym.sklep.pl/system-do-dokumentacji-zeli-z-transiluminatorem-uv-slide-imager-with-uv-transilluminator-p-170.html>



System do dokumentacji żeli z transiluminatorem UV (Slide Imager with UV Transilluminator)

Dostępność

Na zamówienie

Numer katalogowy

SMU-01

Opis produktu

SLIDE IMAGER SMU-01

Features?

- Easy: gel photography and cutting at the same instrument
- Affordable: competitive price
- Quality: branded digital camera,
- ~100 % dark room for photography
- LED trans or UV trans Embedded
- Friendly: USB interface, no need to take out memory card frequently
- Considerate: rotatable LCD monitor for user to easily perform photography
- Safe: cutting gel with MaestroGen UltraSafe UV shield (99.99 % UV cut off)