

Link do produktu: <http://www.novazym.sklep.pl/thermal-cycler-gene-explorer-ge4832-48-x-02-ml-30-x-05-ml-p-333.html>



Thermal Cycler Gene Explorer GE4832 (48 x 0.2 mL + 30 x 0.5 mL)

Dostępność

Na zamówienie

Numer katalogowy

GE-4832

Opis produktu

Thermal Cycler Gene Explorer GE4832 (48 x 0.2 mL + 30 x 0.5 mL)

GENE EXPLORER

The product has a variety of advanced technologies: WINDOWS operating system, color display screen, dual modules operate the higher requirements for an experiment. The quantitative PCR product perfect to the extreme, powerful features will meet

Features:

- 1) The most advanced peltier-based semiconductor technology.
- 2) High-performance universal power supply.
- 3) Large 5.6-inch high-definition color LCD screen.
- 4) Windows operating system, user-friendly interface.
- 5) Graphical user interface in English and Chinese.
- 6) Online features and print capabilities.
- 7) Power-down data protection.
- 8) Friendly design of the shell structure, a strong three-dimensional visual impact.
- 9) Stepless adjustable hot lid.
- 10) Lid can be positioned at any angle.
- 11) High-sealing reaction zone, to ensure stable and reliable test.
- 12) Independent control modules, Can be used at least two at the same time
- 13) Modular body design, easy to maintain and update.
- 14) Eight modules can be replaced easily meet the requirements of different experiments.