

Dane aktualne na dzień: 12-05-2024 17:15

Link do produktu: <http://www.novazym.sklep.pl/thermo-shaker-mix-smart-mix-100-p-136.html>



Thermo-shaker MIX Smart MIX-100

Dostępność

Na zamówienie

Numer katalogowy

MIX-100

Opis produktu

[Thermo-shaker MIX-100 Mix-Smart](#)

{C}{C}



{C}{C}



{C}{C}



{C}{C}

compact and amazingly designed, benchtop Mix-Smart can be used for virtually any mixing applications. It accommodates micro test tubes, PCR plates, deepwell plates and MTPs.

It features an integrated vortexing function and an operating concept that is as simple as it is intelligent. The Mix-Smart only takes up an area of 17 x 22 cm, less than the size of this A4 paper. But the Mix-Smart can be offered many options for holding a wide variety of tubes and plates:

Features:

- Mixing in a 3-in-1 format
- Plates (up to 96- and 384-well formats)
- PCR plates (skirted, semi-skirted, unskirted)
- Deepwell plates, MTPs, Tubes
- 0.2 ml PCR tubes and PCR strips
- 0.5 ml/1.5ml/2.0ml micro test tubes
- Vortexing numerous tube formats
- Preprogrammed soft keys and Reliable universal holder
- Automatic imbalance detection— a common feature in centrifuges
- Integrated vortex function

Specification:

- Power supply: 230V/120V, 50-60 Hz
- Power: 45 W
- Dimensions (WxDxH): 17x22x12 cm
- Weight: 4.5 kg
- Mixing Speed: 300 to 3,000 rpm
- Touch vortexing Speed: 3,500 rpm
- Adjustable mixing time: 15 s to 99.5 h; continuous
- Mixing and vortexing radius: 1.5 mm (3 mm mixing Orbit)