

Link do produktu: <http://www.novazym.sklep.pl/ultrabright-uv-transilluminator-302-nm-p-386.html>



UltraBright UV transilluminator (302 nm)

Dostępność

Na zamówienie

Numer katalogowy

MB-16

Opis produktu

UltraBright UV transilluminator MaestroGen

Specification

- Dimensions (mm): 280D x 340W x 80H
- Viewing surface (mm): 160D x 200W
- Wavelength: 302 nm
- UV tubes: 5 x 8W
- UV blocking cover: Aluminum alloy frame
- Power: 100-240 V, 50-60 Hz
- Weight: 4.3 kg

Key features

- Uniform lighting design
- UltraSafe UV blocking cover: block 99.99% UV transmittance
- UV wavelength: 302nm
- Adjustable intensity: from 100% down to 50%