

Link do produktu: <http://www.novazym.sklep.pl/ultrabright-uv-transilluminator-large-302-nm-large-viewing-p-385.html>



UltraBright UV Transilluminator LARGE (302 nm, large viewing)

Dostępność

Na zamówienie

Numer katalogowy

MB-21

Opis produktu

UltraBright UV Transilluminator MaestroGen

Specification

- Dimensions(mm): 280D x 340W x 80H
- Viewing surface(mm): 210D x 260W
- Wavelength: 302 nm
- UV tubes: 6 x 8 W
- UV blocking cover: Aluminum alloy frame
- Power: 100-240 V, 50-60 Hz
- Weight:: 4.3 kg

Key features

- Uniform lighting design
- UltraSafe UV blocking cover: block 99.99% UV transmittance
- UV wavelength: 302nm
- Adjustable intensity: from 100% down to 50%