

Link do produktu: <http://www.novazym.sklep.pl/ultrabright-uv-transilluminator-large-dual-uv-wavelengths-302-365-nm-large-viewing-p-384.html>



UltraBright UV Transilluminator LARGE (dual UV wavelengths 302 & 365 nm, large viewing)

Dostępność

Na zamówienie

Numer katalogowy

MLB-21

Opis produktu

UltraBright UV Transilluminator LARGE MaestroGen

MaestroGen UltraBright UV transilluminator offers UltraSafe UV shield and strong frame design, no UV shield breakage anymore!

Specification

- Dimensions (mm): 280D x 340W x 80H
- Viewing surface (mm): 210D x 260W
- Wavelength: 302 & 365 nm
- UV tubes: 11 x 8 W
- UV blocking cover: Aluminum alloy frame
- Power: 100-240 V, 50-60 Hz
- Weight : 4.3 kg

Key features

- Uniform lighting design, Ultra High Contrast!
- UltraSafe UV blocking cover: block 99.99% UV transmittance
- Strong frame and stylish design (no breakage of UV shield)
- Fast turn on speed for UV lamps
- Dual UV wavelengths: 302nm and 365nm
- Adjustable intensity: from 100% down to 50%