

Link do produktu: <http://www.novazym.sklep.pl/ultrabright-uv-transilluminator-quadratic-302-nm-quadratic-viewing-surface-p-387.html>



UltraBright UV Transilluminator QUADRATIC (302 nm, quadratic viewing surface)

Dostępność

Na zamówienie

Numer katalogowy

MB-20

Opis produktu

UltraBright UV Transilluminator QUADRATIC MaestroGen

Specification

- Dimensions(mm): 280D x 340W x 80H
- Viewing surface(mm): 200D x 200W
- Wavelength: 302 nm
- UV tubes: 6 x 8 W
- UV blocking cover: Aluminum alloy frame
- Power: 100-240 V, 50-60 Hz
- Weight: 4.3 kg

Key features

- Uniform lighting design
- UltraSafe UV blocking cover: block 99.99% UV transmittance
- UV wavelength: 302nm
- Adjustable intensity: from 100% down to 50%