

Link do produktu: <http://www.novazym.sklep.pl/ultrabright-uv-transiluminator-dual-uv-wavelengths-302-365-nm-p-162.html>



UltraBright UV Transiluminator (dual UV wavelengths 302 & 365 nm)

Dostępność

Na zamówienie

Numer katalogowy

MLB-16

Opis produktu

UltraBright UV Transiluminator MaestroGen

Specification

- Dimensions (mm): 280D x 340W x 80H
- Viewing surface (mm): 160D x 200W
- Wavelength: 302 & 365 nm
- UV tubes: 9 x 8W
- UV blocking cover: Aluminum alloy frame
- Power: 100-240 V, 50-60 Hz
- Weight: 4.3 kg

Key features

- Uniform lighting design, Ultra High Contrast!
- UltraSafe UV blocking: block 99.99% UV transmittance
- Strong and stylish metal frame design (NO breakage of UV shield)
- Fast turn on speed for UV lamps
- Dual UV wavelengths: 302 nm and 365 nm
- Adjustable intensity: from 100% down to 50%