

Link do produktu: <http://www.novazym.sklep.pl/ultraslim-led-transiluminator-black-470-nm-easy-to-cut-gel-p-153.html>



UltraSlim LED Transiluminator BLACK (470 nm, easy to cut gel)

Dostępność

Na zamówienie

Numer katalogowy

SLB-01

Opis produktu

UltraSlim LED Transiluminator BLACK MaestroGen

MaestroGen UltraSlim LED illuminator offers the feasibility of easy gel cutting and stylish design.

Specification

- Dimensions (mm): 210D x 210W x 30H
- Viewing surface (mm): 120D x 70W
- Sensitivity: up to 0.5 ng
- Uniformity: < 10 % CV (coefficient of variance)
- LED Wavelength: 470 nm
- Power: 24 V - 1.25 A
- Weight: 2.3 kg

Key features

- Compact: 30 mm in height, THINNEST in the world now to enlarge your laboratorial space, 210(D)X210(W)X30(H) mm
- User Friendly: Easy to view and cut gel (NO need to wear amber glasses)
- Light: Convenient to move around (for personal)
- Safe: Blue LEDs do not cause burn to your eyes and skins; NO DNA damage inherent to UV
- Broad application : Optimized for use with most nucleic acid and protein stains
- Uniform: supplies <10% coefficient of variance (CV) in the central viewing area