

Dane aktualne na dzień: 02-05-2025 05:32

Link do produktu: <http://www.novazym.sklep.pl/uv-ii-uv-portable-transilluminator-p-372.html>



UV-II UV Portable Transilluminator

Dostępność

Na zamówienie

Opis produktu

UV-II UV Portable Transilluminator

JY series UV Transilluminator is mainly applied to the electrophoresis results observation. With the principles of different wavelength shows the different electrophoresis band result, They are compatible with the observation of nuclear acid and protein electrophoresis bands, by using the high resolution of CCD and software technology, the gel pictures are clearly showed.

Specification

- Reflection ultraviolet wavelength: 254 nm 365 nm
- Two options for alternate use