

Link do produktu: <http://www.novazym.sklep.pl/vetpcr-p-haemolytica-48d-p-740.html>



VetPCR? P. haemolytica (48D)

Numer katalogowy

VET-B020-48D

Opis produktu

Pasteurella mannheimia haemolytica PCR Detection Kit

Mannheimia haemolytica is the principal bacterium isolated from respiratory disease in feedlot cattle and is a significant component of enzootic pneumonia in all neonatal calves. A commensal of the nasopharynx, M. haemolytica is an opportunist, gaining access to the lungs when host defenses are compromised by stress or infection with respiratory viruses or mycoplasma. M. haemolytica may contribute to enzootic pneumonia of calves; enzootic pneumonia of lambs and peritonitis in sheep. It also causes gangrenous mastitis in sheep.

VetPCR? M. haemolytica Detection Kit is the direct detection of Pasteurella Mannheimia haemolytica on the basis of a genetic database, so it can diagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 3 hours for detection. Therefore, it is a very fast, accurate and reliable technique.

Characteristics

Ready to use : only DNA template and D.W. are needed. Easy and speed protocol. Stable for 1 year at -20°C. Time-saving and cost-effective.

Contents

Kit	Quantity 48	Quantity 96	Package
VetPCR™ M. haemolytica Premixture	1	1	Vial
PCR Internal Control	1	1	Vial
DNase/RNase free water	1	1	Vial
M. haemolytica PCR Positive control	1	1	Vial
PCR Negative control	1	1	Vial
Mineral Oil Solution	1	1	Vial
Brig™ Molecular Weight marker	1	1	Vial
DNA purification kit	50	100	Test

Interpretation of the Test Result

Expected PCR product size : 259bp



Fig. 1 Result:

Lane M: Brig? Molecular Weight Marker

Lane 1~2: M. haemolytica Positive samples

Lane I.C.: Internal control

Lane P: Positive control

Lane N: Negative control