

Link do produktu: <http://www.novazym.sklep.pl/biological-indicator-incubator-minit-1-p-518.html>



# Biological Indicator Incubator MiniT-1

Dostępność

Dostępny

## Opis produktu

### Biological Indicator Incubator Model? MiniT-1

The MiniT-1 Biological Indicator Incubator is designed for temperature incubation of the EO Biological Indicators. It provides easy verification of incubation temperature.

#### Features

- 1. LCD display, with 12V power for ease of portability
- Operating Temperature Range: From 5°C above ambient to 80°C
- Time Range: 1min ~ 99h59min
- 4Temp.Control Accuracy:  $\leq \pm 0.5^\circ\text{C}$
- 5Display Accuracy:  $0.1^\circ\text{C}$
- 6Heating Time: < 6min (from 20°C to 57°C)
- Temp.Control Accuracy @57°C:  $\leq \pm 0.35^\circ\text{C}$
- Dimensions: 110(w) x 150(D) x 80(H) mm

#### Specification

- 1. The Incubator has many wells for placement of the biological indicator vials and the wells can be custom-made. User can set the incubating temperature.
- 2. Materials: Aluminium block
- 3. Well Diameter: According to custom-making

| Type      | Description       | Remark |
|-----------|-------------------|--------|
| A-MiniT-1 | 15 x $\phi$ 10mm  |        |
| B-MiniT-1 | 12 x $\phi$ 12mm  |        |
| C-MiniT-1 | 6 x $\phi$ 15mm   |        |
| D-MiniT-1 | 6 x $\phi$ 16mm   |        |
| E-MiniT-1 | 18 x $\phi$ 8.8mm |        |