

Dane aktualne na dzień: 05-06-2024 02:53

Link do produktu: <http://www.novazym.sklep.pl/mini-clinical-centrifuge-hwlab-p-3811.html>



Mini clinical centrifuge HWLAB

Opis produktu

Mini clinical centrifuge HWLAB



Description

HWLAB® mini clinical centrifuge is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Technical specifications

Description	HW6M
Maximum rcf	2500xg
Power supply	AC100V-240V/50-60Hz
Rotor capacity	6 x 15ml/10ml/7ml/5ml tubes
Speed	From 300r/min to 5000r/min
Timer	15s-99min / continuous operation
Motor type	Brushless motor
Dimension	212x270x158mm
Weight	4.5KG

Product features

- Maintenance-free brushless motor
- LCD display of time and speed
- Adjustable speed from 300r/min to 5000r/min
- Timer, from 15s to 99min, or continuous operation
- For 6x15ml/10ml/7ml/5ml tubes
- Lid automatically open at the end of the run to prevent sample warming and to reduce handling time



Ordering information: HWLAB® Mini clinical centrifuge

	Model	Cat.NO	Destription
	HW6M	C-12000	Widely ues in clinical labs and donation centers
	HW6M.1	C-12001	For 6x15ml/10ml/7ml/5ml tubes