

Link do produktu: <http://www.novazym.sklep.pl/miniwirowka-laboratoryjna-6kc-mini-centrifuge-p-458.html>



## Miniwirówka laboratoryjna 6KC (Mini-Centrifuge)

Dostępność

**Dostępny**

Numer katalogowy

**MC-6KC**

### Opis produktu

Mini-6KC Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a "personal" mini centrifuge. It includes 2 quick release interchangeable rotors, two speed options (4,000rpm & 6000rpm). The electronic brake can reach quick deceleration to reduce handling time.

#### Features

1. New rotor design for increased sample capacity
2. Quick release rotor hub allows for easy rotor exchange
3. The electronic braking system for rapid sample processing
4. Two speed options: 4,000 & 6,000rpm

#### Specification

1. Relative centrifugal force: 850 x g (4,000rpm) / 2,000 x g (6,000rpm)
2. Sample capacity: 8x1.5/2.0 ml tube rotor & 4 place PCR Strip rotor
3. Speed of rotation: 4,000rpm & 6,000rpm
4. Maximum noise: ≤56Db
5. Dimensions: 135 x 147 x 107mm (D x W x H)

#### Two speed options:

BRAK  
ZDJĘCIA



#### Two rotors:

BRAK  
ZDJĘCIA



| Code        | Description                     |
|-------------|---------------------------------|
| AS-08020-00 | Mini-6KC AC100V~240V 50~60Hz 1A |
| AS-08021-01 | 8-1.5/2.0ml angle rotor         |
| AS-08021-02 | 4 -8x0.2ml strip rotor          |
| AS-08061-05 | Adapters for 0.5~0.6ml tubes    |
| AS-08061-06 | Adapters for 0.2ml tubes        |