

Link do produktu: <http://www.novazym.sklep.pl/nextractor-nx-48-p-4134.html>



Nextractor (NX-48)

Numer katalogowy

NX-48

Opis produktu

Nextractor (NX-48)

Nextractor is automatic extraction system for rapid isolation of DNA or RNA from a variety samples.

Using of Nextractor can purify up to 48 samples, saving valuable time and labor in your lab.

Also it offers a high throughput processing and yields.

Specification:

1. Fully Automated Extraction

-Walk away system

2. High Throughput

-48 Specimens at one run

-48 Extractions in 10 minutes

3. Pre-filled Catridge type of Reagent

-Single, 8, 24 Extraction Cartridge

-All Cartridge can be used simultaneously without and change of Hardware and software format

4. Safety

-UV Lamp/Auto door Lock

5. Touch Screen based User interface

-Easy to operate without training

-Pre-programmed and User define protocols

6. Dimension

-364mm(W)x386mm(H)x420mm(D)

-25kg