

Link do produktu: <http://www.novazym.sklep.pl/wytrzasarka-os-300-orbital-shaker-p-298.html>



Wytrząsarka OS-300 (Orbital Shaker)

Dostępność

Na zamówienie

Numer katalogowy

OS-300

Opis produktu

OS-300 Orbital Shaker

The OS-300 Orbital Shaker holds a Multi-function Shaking system. And the intensive multi-function shaking system surpasses all previous levels of experience! You can set the shaking modes in maximum nine ways of Forward / Backward / Pause / Mix & match upon your various laboratory purposes.

Specification

Shaking speed	50 – 250 rpm
Orbit	10, 15, 20 mm available (15 mm as Factory default)
Independent timer	0 – 99 h 59 min
Load capacity	8 kg
Net weight	17.5 kg
Dimensions (mm)	320 x 350 x 145
Base platform size (mm)	270 x 240
Program function	Programmable, 9 steps, 9 programs
Platform	Universal Platform (see below)
Flask clamps	see below
Motor	Brushless motor
Power supply	AC80 V ~ 240 V, 50/60 Hz

Key features

- Professional control electronics and programmable control (microprocessor PID controls, programmable control; self-diagnostic function identifies any errors)
- Multi-function shaking system (run time; forward & Backward, Pause)
- Brushless DC motor applications for smooth shaking
- Wide range of accessories
- Multi-function Shaking System - intensive multi-function shaking system surpasses all previous levels of experience! You are able to easy set the shaking modes in maximum nine ways of Forward /Backward / Pause / Mix & match upon your various laboratory purpose