

Link do produktu: <http://www.novazym.sklep.pl/bandpeeper-for-mupid-2plus-series-real-time-electrophoresis-p-160.html>



BandPeeper for Mupid 2plus series (Real-Time Electrophoresis)

Dostępność

Na zamówienie

Numer katalogowy

RT-160

Opis produktu

[BandPeeper for Mupid 2plus series \(Real-Time Electrophoresis\) MaestroGen](#)

[Specification](#)

Size: 166*170*82 mm
Window size: 150*60 mm
Filter type: Amber filter
Filter surface: Anti-Fog treatment
PCBA surface: Anti-moisture treatment
Power: 24V - 1.25A
Sensitivity: ~ 2 ng
LED wavelength: 470 nm
Compatibility: only with Mupid 2plus tanks
Weight: 0.96 kg

[Key features](#)

- Real-Time Fluorescent signal observation
- Strong blue LED: used with common nucleic acid dyes
- High sensitivity
- Anti-fog design on amber filter
- Time saving for trial & error of new sample electrophoresis
- Accurate for electrophoresis duration control
- Strongly suggested to bundle with electrophoresis tank, welcome ODM/OEM for special tank